

EARLY EXPERIENCE WITH PATIENT TREATMENT AND DRUG DELIVERY WITH IRRAflow

An Automatically Irrigating and Draining Ventricular Catheter

Presented by Babar Khan, MD

and from Nicholas Brandmeir, MD,
and WVU Rockefeller Neuroscience Institute
at The Neurosurgical Society of the Virginias Annual Meeting

The need for improving EVDs

EVD COMPLICATION RATES

- Infection: 4-20%
 - Increases with occlusion
 - Current abx impregnated catheters only prevent gram + infection
- Revision: 13%
 - Occlusion: 25-50% in IVH
- Hemorrhage: 8%

Conclusion

Preliminary data confirms that IRRAflow might reduce rates of catheter occlusion, infections, symptomatic and radiographic vasospasm.

Future prospective (randomized) studies will be conducted to confirm these findings.

Yuen J, Selbi W, Muquit S, Berei T. Complication rates of external ventricular drain insertion by surgeons of different experience. *Ann R Coll Surg Engl* [Internet] 2018;100:221-5

Atkinson RA, Fikrey L, Vail A, Patel HC. Silver-impregnated external-ventricular-drain-related cerebrospinal fluid infections: a meta-analysis. *J Hosp Infect* [Internet] 2016;92:263-72

Talibi SS, Silva AH, Afshari FT, Hodson J, Roberts SA, Oppenheim B, et al. The implementation of an external ventricular drain care bundle to reduce infection rates. *Br J Neurosurg* [Internet] 2020;34:181-6.

Gilard V, Djoubairou BO, Lepetit A, Metayer T, Gakuba C, Gourio C, et al. Small versus Large Catheters for Ventriculostomy in the Management of Intraventricular Hemorrhage. *World Neurosurg* [Internet] 2017;97:117-22.



Applications

- **IVH**
 - Irrigation
 - tPA infusion
- **VENTRICULITIS**
 - IT infusion of abx
- **SAH**
 - IT infusion of anti-vasospasm agents (nicardipine)
- **SDH**
 - Irrigation




Patient Characteristics

45 patients

Gender	Male	28 (62%)
	Female	17 (38%)
Length of Stay (days)	ICU	7.4 ± 5.8
	Hospital	14.4 ± 9.8
	Age (years)	60.2 ± 17.3
Indications	IVH	23 (51.1%)
	SAH HH1 (0), HH2 (2), HH3 (6), HH4 (2), HH5 (2)	12 (26.6%)
	SDH	13 (28.8%)
	Ventriculitis	3 (6.6%)
	Hydrocephalus (Other)	2 (4.4%)
Insertion Location	R frontal	26 (57.7%)
	L frontal	6 (13.3%)
	Burr hole	5 (11.1%)
	Craniotomy	8 (17.7%)
Average Duration	(days)	6.8 ± 4.9
Insertion Method	Stealth (bedside)	18 (40%)
	Anatomic	14 (31.1%)
	Craniotomy/burrhole	13 (28.8%)


CT and endoscopic images confirm impact of IRRAflow's automated irrigation



Pre-IRRAflow

Post-IRRAflow

- IRRAflow inserted directly in clot
- Irrigation performed with tPA
- Endoscopic evaluation performed on day 4



NO OCCLUSION FORMATION EVEN IN HEAVY CLOT BURDEN

Irrigation from IRRAflow cleansing the catheter tip, keeping the **catheter free of occlusions**

Images provided by Dr. Nicholas Brandmeir

Outcomes

Complications	Ventriculitis	1 (2.2%)
	Mortality	9 (20%)
Intrathecal Medications administered	tPA	9 (20%)
	Vancomycin	2 (4.4%)
Catheter Occlusion		0%
Shunt dependence (IVH)		3/23 (13%)**
Vasospasm	Clinical	2/12 (17%)***
	Radiographic	3/12 (25%)***
SDH	MLS before: 4.9±2.6	Percentage improved 13/13 (100%)
	MLS after 2d: 0.8±1.2	

* Compared to 19% permanent, 41% temporary occlusion.

Fargen et al. JNS 2016

** Compared to 18% in literature. Clear III trial.

Murthy et al. Neurology 2017

*** Compared to 40-70% in literature.

Bracard et al. Interventional Neuroradiology 2008